



Inspect these 3 points monthly.



Peel here (Recommended locations: nightstand, bedside table, bathroom mirror or refrigerator)

Peel here (Adhere to your CPAP machine)



Use these guidelines to help you inspect, assess and replace your CPAP supplies. Ask your supplier how.

EVERY 1 MONTH

Mask cushions and pillows

- Have the edges of your cushion or pillows become stiff or cracked?
- Has your cushion changed color?
- Are you experiencing discomfort, redness or skin irritation?

Machine filters

- Does your filter show signs of wear such as discoloration, holes or tears?

EVERY 3 MONTHS

Mask (excludes headgear)

- Are you experiencing excessive leaks?

Tubing (standard or heated)

- Has your tubing begun to change to an opaque color?
- Has your tubing developed tears?

EVERY 6 MONTHS

Headgear and chin strap

- Has your headgear stretched?
- Do you find that you are over-tightening to avoid leaks?

Humidifier water chamber

- Has the water chamber become discolored or cloudy?
- Do you see any cracks or pitted areas?

Visit ResMed.com/register to stay on track with your therapy.

© 2015 ResMed. Specifications may change without notice. 1016060/3 2015-04



ResMed