Device setting
Full face

AirTouch™
F30i
COMFORT

AirTouch™
F30i
CLEAR


AirFit™
F30i

User Guide

 **Resmed**


English | Bahasa Melayu
Indonesia | 繁體中文 | 日本語
한국어 | ไทย | Tiếng Việt

- | | |
|----------------------|--|
| A Elbow | 5 Anti-Asphyxia valve (inside the left and right of frame connector) |
| B Frame | 6 Magnetic clip |
| C Headgear | 7 Lower headgear strap |
| 1 Vent holes (Elbow) | 8 Upper headgear strap |
| 2 Side button | 9 Frame magnet |
| 3 Elbow ring | |
| 4 Frame connector | |

AirTouch™ F30i


Full face mask

- D Cushion
- 10 QuietAir™ vent
- 11 Multi-hole vent

AirFit™ F30i

Full face mask

- E Cushion
- 12 Vent holes (Cushion)

Notes:

- The AirTouch F30i full face mask features a fabric cushion and is available in two frame options:
 - The AirTouch F30i Comfort features a fabric-wrapped frame. The fabric sleeve is not removable.
 - The AirTouch F30i Clear features a silicone frame.
- The AirFit F30i full face mask features a silicone cushion and a silicone frame.

- The cushions for these masks are interchangeable.
- Not all masks are available in all regions.

Intended use

The AirTouch F30i and AirFit F30i masks are intended to be used by patients weighing more than 30 kg who have been prescribed non-invasive positive airway pressure (PAP) therapy such as CPAP or bi-level therapy. The mask is intended for single patient re-use in the home and multi-patient re-use in the hospital/institutional environment.

Contraindications

Masks with magnetic components are contraindicated for use by patients where they, or anyone in close physical contact while using the mask, have the following:

- Active medical implants that interact with magnets (ie, pacemakers, implantable cardioverter defibrillators (ICD), neurostimulators, cerebrospinal fluid (CSF) shunts, insulin/infusion pumps)
- Metallic implants/objects containing ferromagnetic material (ie, aneurysm clips/flow disruption devices, embolic coils, stents, valves, electrodes, implants to restore hearing or balance with implanted magnets, ocular implants, metallic splinters in the eye).

WARNING

Keep the mask magnets at a safe distance of at least 150 mm away from implants or medical devices that may be adversely affected by magnetic interference. This warning applies to you or anyone in close physical contact with your mask. The magnets are in the frame and lower headgear clips, with a magnetic field strength of up to 400 mT. When worn, they connect to secure the mask but may inadvertently detach while asleep.

Implants/medical devices, including those listed within contraindications, may be adversely affected if they change function under external magnetic fields or contain ferromagnetic materials that attract/repel to magnetic fields (some metallic implants, eg, contact lenses with metal, dental implants, metallic cranial plates, screws, burr hole covers, and bone substitute devices). Consult your physician and manufacturer of your implant/other medical device for information on the potential adverse effects of magnetic fields.

WARNING

- The mask contains safety features, the exhaust vent holes and anti-asphyxia valves, to enable normal breathing and exhaust exhaled breath. Occlusion of the exhaust vent holes or anti-asphyxia valves needs to be prevented to avoid having an adverse effect on the safety and quality of the therapy. Regularly inspect the vent holes and anti-asphyxia valves to ensure they are kept clean, clear of blockages and are not damaged.
- If there is any visible deterioration of a mask component (cracking, crazing, tears etc.) the component should be discarded and replaced.
- Only use compatible CPAP or bi-level therapy devices or accessories. The technical specifications of the mask are provided for healthcare professionals to determine compatible devices. Use in combination with incompatible medical devices can decrease the safety or alter the performance of the mask.
- The mask is not suitable for patients requiring life support ventilation or who would experience serious deterioration of health or death with a loss or degradation of therapy.
- Regularly clean your mask and its components to maintain the quality of your mask and to prevent the growth of germs that can adversely affect your health.
- The mask must be used under qualified supervision for patients who are unable to remove the mask by themselves. The mask may not be suitable for those prone to aspiration.
- The mask should not be used unless the device is turned on. Once the mask is fitted, ensure the device is blowing air to reduce risk of rebreathing exhaled air.
- Discontinue using or replace this mask if the patient has ANY adverse reaction to the use of the mask. Consult your physician or sleep therapist.
- Follow all precautions when using supplemental oxygen.
- Oxygen flow must be turned off when the CPAP or bilevel device is not operating, so that unused oxygen does not accumulate within the device enclosure and create a risk of fire.
- Oxygen supports combustion. Oxygen must not be used while smoking or in the presence of an open flame. Only use oxygen in well ventilated rooms.
- At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration varies, depending on the pressure settings, patient breathing pattern, mask, point of application and leak rate. This warning applies to most types of CPAP or bilevel devices.

WARNING

- The mask is not intended to be used simultaneously with nebulizer medications that are in the air path of the mask/tube.
- The mask is Magnetic Resonance (MR) unsafe and must be kept outside of MRI scanner rooms.
- Always follow cleaning instructions and only use a mild liquid detergent. Some cleaning products may damage the mask, its parts and their function, or leave harmful residual vapours. Do not use a dishwasher or washing machine to clean the mask. Ozone or UV light products have not been validated for use with the mask and may lead to discoloration, damage or reduced mask seal performance.

CAUTION


- When fitting the mask, do not overtighten the headgear as this may lead to skin redness or sores around the mask cushion.
- As with all masks, some rebreathing may occur at low pressures.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult your physician or dentist.
- Do not iron the headgear as the material is heat sensitive and will be damaged.

Note: For any serious incidents that occur in relation to this device, these should be reported to Resmed and the competent authority in your country.

Before using your mask

Remove all packaging before using the mask. Wrinkles on fabric components may occur during storage and transportation, and do not affect the performance or safety of the mask.

Fitting your mask

1. Connect the air tubing from your device to the elbow.
2. Squeeze the side buttons on the elbow and pull from the frame. Put the elbow and air tubing aside for now.
3. Twist and pull both magnetic clips away from the frame magnets.
4. Place the cushion under your nose and ensure it sits comfortably against your face. With the Resmed logo on the headgear facing up, pull the headgear and frame over your head.
5. Bring the lower headgear straps under your ears and attach the magnetic clips to the frame magnets.
6. Undo the fastening tabs on the upper headgear straps and pull evenly. Repeat with the lower headgear straps.

7. Attach the elbow to the top of the frame. Your mask should now be positioned as shown.

Adjusting your mask

- When using the mask, air will flow out of the vents in the cushion and elbow. If air leaks from the side or top of the cushion, adjust the mask position to improve the seal.
- To fix any mask leaks, adjust the upper or lower headgear straps. Adjust only enough for a comfortable seal and do not overtighten.

Scan this QR code with your mobile device or visit [Resmed.com/downloads/masks](https://resmed.com/downloads/masks) for helpful documents and support videos for your mask.

Mask setting options

Ensure that **Full Face** is selected in the mask setting options of your Resmed device. For a full list of compatible devices for this mask, see the Mask/Device Compatibility List at [Resmed.com/downloads/masks](https://resmed.com/downloads/masks).

Device Setting

Full Face


Sizing the mask cushion

If available, your physician or sleep therapist will use the fitting template to assist in selecting the right size cushion based on nose shape. Note that sizing across different masks may not be the same.

	AirTouch F30i	AirFit F30i
Narrower noses	M	S or M
Wider noses	SW or L	SW or W

Sizing the mask frame

If the mask frame falls back over your head or sits too close to your ears, try a smaller frame size. If the mask falls forward on your head or sits too close to your eyes, try a bigger frame size.


Optimizing mask seal and comfort


To receive the most effective quality of therapy it is important to have the best seal between your mask and your face. If you are experiencing leaks ensure that you have followed the fitting instructions carefully. After fitting your mask check the following:

- The headgear and frame is fitted correctly. Readjust the tension of the headgear and reseat the frame position to ensure a comfortable fit and proper seal.
- The correct cushion sizes are used. Try a larger or smaller cushion.
- The mask is clean and free of oils on the cushions.

Disassembling your mask for cleaning

If your mask is connected to a device, disconnect the device air tubing from the mask elbow.

1. Undo the fastening tabs on the upper headgear straps and pull from the frame. Keep the magnetic clips attached to the lower headgear straps.
2. Squeeze the side buttons on the elbow and detach from the frame.
3. Grip the frame connector and lift to unclip from the cushion. Repeat on the other side.

Note: After detaching the headgear, close the fastening tabs to prevent it from damaging the fabric sleeves when handling.

Cleaning your mask at home

If there is any visible deterioration of a mask component (cracking, crazing, tears etc.) the component should be discarded and replaced. Pilling and fraying of fabric components may occur and do not affect performance of the mask.

Thoroughly clean the mask components as follows


Daily/After each use: Cushion

Weekly: Headgear, Frame, Elbow

1. Soak the components in warm water with a mild liquid detergent.
2. Hand wash the components with a soft, non-metallic bristle brush (eg, soft bristle toothbrush). Pay particular attention to the vents in both cushion and elbow.
For the AirTouch F30i cushion, hand wash only without a bristle brush.
3. Thoroughly rinse the components under running water.
For AirTouch F30i Comfort, if the sleeve on the frame becomes twisted, it can be repositioned when wet.
4. Squeeze the fabric components with a clean towel to remove excess water.
5. For the QuietAir vent, shake to remove any excess water.
6. Leave the components to air dry out of direct sunlight.

If the mask components are not visibly clean, repeat the cleaning steps. Make sure that the vents and anti-asphyxia valves are clean and clear.

Reassembling your mask


1. Align and insert the frame connector tab (a) into the cushion slot (b) and press down until it clicks. Repeat on the other side.
2. Attach the elbow to the top of the frame until it clicks.
3. With the Resmed logo facing outwards and up, insert the upper headgear straps into the frame from the inside
4. Fold the fastening tabs over to secure.


Notes:

- If the elbow ring detaches, re-insert into the top of the frame.
- For AirTouch F30i, if the vent detaches, refer to the **Removing and inserting the vent** section of this guide.

Removing and inserting the vent (AirTouch F30i)

1. To remove the vent, gently peel the cushion away and slide the vent out.

2. To insert the vent, align the tab (a) on the cushion with the notch (b) on the vent, and then gently pull the cushion over the vent.

Reprocessing the mask between patients


When using between patients, these masks must be reprocessed according to the cleaning and disinfection instructions available on [Resmed.com/downloads/masks](https://www.resmed.com/downloads/masks).

Technical specifications

Mask setting options: For AirSense™, AirCurve™ or S9™ devices: Select 'Full Face'.

Compatible devices: For a full list of compatible devices for this mask, see the Mask/Device Compatibility List at [Resmed.com/downloads/masks](https://resmed.com/downloads/masks). SmartStart/SmartStop may not operate effectively when using this mask with some CPAP or bi-level devices.

Pressure-flow curve

Mask pressure (cmH ₂ O)	AirTouch F30i Vent flow rate (L/min)	AirFit F30i
4	22	22
11	39	38
17	49	48
24	60	59
30	69	67
Therapy pressure	4 - 30 cmH ₂ O	4 - 30 cmH ₂ O
Resistance with Anti Asphyxia Valve closed to atmosphere		
Drop in pressure measured (nominal)		
at 50 L/min	0.3 cmH ₂ O	0.2 cmH ₂ O
at 100 L/min	1.2 cmH ₂ O	1.0 cmH ₂ O
Resistance may vary due to the flexible mask frame design.		
Resistance with Anti Asphyxia Valve open to atmosphere		
Inspiration at 50 L/min	0.2 cmH ₂ O	0.2 cmH ₂ O
Expiration at 50 L/min	0.3 cmH ₂ O	0.4 cmH ₂ O
Anti-Asphyxia Valve		
open-to-atmosphere pressure	<4 cmH ₂ O	<4 cmH ₂ O
closed-to-atmosphere pressure	<4 cmH ₂ O	<4 cmH ₂ O

	AirTouch F30i	AirFit F30i
Sound		
Declared dual-number noise emission values in accordance with ISO4871:1996 and ISO3744:2010. A-weighted sound power level and A-weighted sound pressure level at a distance of 1 m, with uncertainty of 3 dBA, are shown		
Sound power level	-	25 dBA
Sound power level (with QuietAir vent)	20 dBA	-
Sound power level (with Multi-hole vent)	24 dBA	-
Sound pressure level	-	18 dBA
Sound pressure level (with QuietAir vent)	13 dBA	-
Sound pressure level (with Multi-hole vent)	17 dBA	-
Gross dimensions (H x W x D)		
Mask fully assembled, excluding headgear (+/- 10% on all measurements)	298 mm x 162 mm x 67 mm	307 mm x 161 mm x 62 mm

Environmental conditions

Operating temperature: 5°C to 40°C

Operating humidity: 15% to 95% RH non-condensing

Storage and transport temperature: -20°C to +60°C

Storage and transport humidity: up to 95% RH non-condensing

International Commission on Non-Ionizing Radiation Protection (ICNIRP)

Magnets used in this mask are within ICNIRP guidelines for general public use. The static magnetic field strength is less than 400 mT at component surface and less than 0.5 mT at 50 mm distance.

Service life: The service life of the mask system is dependent on the intensity of usage, maintenance, and environmental conditions to which the mask is used or stored. As this mask system and its components are modular in nature, it is recommended that the user maintain and inspect it on a regular basis, and replace the mask system or any components if deemed necessary or according to the instructions in the 'Cleaning your mask' section of this guide.

Storage

Ensure that the mask is thoroughly clean and dry before storing it for any length of time. Store the mask in a dry place out of direct sunlight.

Disposal

This mask and packaging does not contain any hazardous substances and may be disposed of with your normal household refuse.

Symbols

The following symbols may appear on your product or packaging.

Full face mask


Device setting - Full face

Small frame

Standard frame

See symbols glossary at [Resmed.com/symbols](https://resmed.com/symbols).

Limited warranty

Our goods come with guarantees that cannot be excluded under the Australian and New Zealand Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

In addition to your rights and remedies under the Australian and New Zealand Consumer Law (and any other applicable law), ResMed Pty Ltd ABN 30 003 765 142 of 1 Elizabeth Macarthur Drive, Bella Vista NSW 2153, (Resmed) warrants that your Resmed product will be free from defects in material and workmanship from the date of purchase for the period specified below:

Product	Warranty period
- AirTouch™ UltraSoft™ memory foam cushion - AcuCare™ series	7 days
- Mask systems (including mask frame, cushion, headgear and elbow) – excluding Consumables - Accessories – excluding Consumables - Flex-type finger pulse sensors - Humidifier water tubs (non-reusable)	90 days
Batteries for use in Resmed internal and external battery systems	6 months
- Clip-type finger pulse sensors - CPAP and bilevel device data modules - DC/DC Converters - Oximeters and CPAP and bilevel device oximeter adapters - Humidifiers and humidifier water tubs (reusable) - Titration control devices	1 year

Product	Warranty period
• CPAP, bilevel and ventilation devices (including external power supply units and excluding humidifier tubs) • Battery accessories (including but not limited to DC cable, PSU adapter and coupler kit) • Portable diagnostic/screening devices	2 years

To make a claim under this warranty you should contact the Resmed Authorised Dealer from which you purchased your Resmed product or contact sleepvantage via phone 1300 305 705 (AUS), 0800 737 633 (NZ) or email: info@sleepvantage.com.au. All claims under this warranty must be accompanied by your original receipt.

You will then need to deliver the Resmed product you claim is defective to the Resmed Authorised Dealer from which you purchased your Resmed product or your closest Resmed Authorised Dealer at your expense. A similar product will normally be lent to you by your Resmed Authorised Dealer while your product is assessed.

The product you claim as defective must be delivered from the Resmed Authorised Dealer to the Resmed Service Centre within the relevant warranty period referred to above. Resmed will not be responsible for the cost of the transport of your Resmed product to the Resmed Service Centre. You must pay any necessary costs to the Resmed Authorised Dealer. If Resmed determines that your warranty claim is valid, we will return the repaired product, or a replacement product, to your Resmed Authorised Dealer at Resmed's expense. If Resmed determines that your warranty claim is valid you may claim any reasonable expenses you have incurred in making the claim by contacting sleepvantage via phone 1300 305 705 (AUS), 0800 737 633 (NZ), email: info@sleepvantage.com.au or posting to us at 1 Elizabeth Macarthur Drive, Bella Vista, NSW 2153 a claim in writing attaching original receipts for the expenses claimed. If Resmed determines that your warranty claim is not valid, we will notify your Resmed Authorised Dealer by providing a quotation of the cost of repair which may include cost of transport of returning your Resmed product back to your Resmed Authorised Dealer. Your Resmed Authorised Dealer will then contact you and you will have the option of taking up the quotation offer, valid for 30 days, or have your product returned unrepared to your Resmed Authorised Dealer at Resmed's expense.

This limited warranty is only available to the initial consumer. It is not transferable.

This limited warranty does not cover: a) any damage caused as a result of improper use, abuse, modification or alteration of the product; b) repairs carried out by any service organisation that has not been expressly authorised by Resmed to perform such repairs; c) any damage or contamination due to cigarette, pipe, cigar or other

smoke and d) any damage caused by exposure to ozone, activated oxygen or other gasses.

This limited warranty is void on product sold, or resold, outside the region of original purchase. Limited warranty claims on defective product must be made by the initial consumer at the point of purchase or through sleepvantage as specified above.


Resmed shall not be responsible for any incidental or consequential damages claimed to have resulted from the sale, installation or use of any Resmed product.

This warranty gives you specific legalrights. For further information on your warranty rights, contact your local Resmed Authorised Dealer or Resmed office.

If you are provided with a replacement or repaired product, the warranty continues to apply to the replacement or repaired device but does not continue beyond the original warranty period referred to on the previous page.

Visit [Resmed.com](https://www.resmed.com) for the latest information on Resmed's Limited Warranty.

If you have any questions or would like the address of your nearest Resmed Authorised Dealer, please contact our friendly sleepvantage team.

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© 2025 ResMed. 638554/1 2025-07
See Resmed.com for other Resmed locations worldwide. AirFit, AirTouch, ComfiSoft, UltraSoft, QuietAir, AcuCare, AirSense, AirCurve, and S9 are trademarks and/or registered trademarks of the Resmed family of companies. For patent and other intellectual property information, see Resmed.com/ip.